

NucleoZOL Total RNA Isolation Protocol


(Xinyi Yu, commented by TCH @ 08/21/2017)

NucleoZOL Reagent, Cat# 740404.200 (Takara BIO, USA)

- 1) **LYSE CELLS:** For cells grown in monolayer, remove cell culture medium and lyse cells by addition of **1mL of NucleoZOL** to the culture dish (e.g., **1ml per one confluent 100mm dish; or 500µl per 60mm dish**). For cells grown in suspension, spin down cells and lyse directly by the addition of NucleoZOL at **1mL NucleoZOL/10⁷ cells**. Lyse cells by pipetting up and down several times. Ensure complete lysis by repeated pipetting.
- 2) **PRECIPITATE PROTEINS/CELL DEBRIS:** Add **200µL** RNase-free water per **500uL NucleoZOL** to the lysate. Shake the sample vigorously for 15 seconds. Centrifuge samples for 5~8min at 12,000 x g at 4°C or room temperature.
- 3) **PC-8 EXTRACTION:** Transfer the RNA-containing supernatant to 2.0ml-Eppendorf tubes, add equal volume of PC-8, vortex vigorously, and centrifuge the mix at top speed for 3-5min. If visible precipitates appear at the interphase, repeat the PC-8 extraction till the interphase is clear. **[NOTE: while this step is optional, it is strongly recommended to avoid any significant RNA degradation, especially when RNA is isolated from tissue samples. PC-8 extraction will remove RNase completely].**
- 4) **PRECIPITATE TOTAL RNA:** Carefully pipette RNA containing top phase from **Step 3** into a 1.7ml or 2.0ml RNase-free Eppendorf tube(s). Add **1mL isopropanol** per **1mL supernatant (i.e., equal volume)** in order to precipitate RNA. **Make sure tubes are kept on ice**. Centrifuge samples for **5~8min** at **12,000 x g** at 4°C or room temperature. **(NOTE: You can stop here by storing the isopropanol-RNA mix at -80°C for days to years).**
- 5) **WASH RNA:** Use **400–600µL 70 % ethanol** when precipitating in 1.5mL or 2mL tubes. For larger tubes, add **500-700uL 70 % ethanol per 1mL supernatant**. Centrifuge the pellets for 1~3 min at 12,000 x g.
- 6) Remove ethanol from the pellet by **carefully pipetting**. Repeat the ethanol washing step. Do not dry the RNA pellet between wash steps.
- 7) **RECONSTITUTE RNA:** Dissolve the RNA pellet in RNase-free water to obtain an RNA concentration of 1~2µg/µL (**i.e., usually dissolve in 20~30µl for one 100mm dish, or 10~15ul for one 60mm dish**). Vortex the sample for 2~5min at room temperature to dissolve total RNA. **Use 1~2µl RNA sample to run RNA “bleach gels” to check the quality/integrity of the sample** (see **Protocol C12** and **see the images below**). Typically, use **5~10µg RNA** per RT-PCR reaction.
- 8) Your RNA samples should be aliquoted and kept at -80°C (**not your own -20°C freezer space**).

NOTE: Similar to TRIzol reagent, NucleoZOL reagent contains a very low level of quinidine to inhibit RNase activity. Thus, you should take the following considerations in mind:

- 1). You should try to collect one sample at a time (keep the rest in the incubator). Once you lyse the cells, please collect the lysate and transfer to 2ml-Eppendorf tubes. Please keep the tubes on ice.
- 2). The RNA will severely degrade if the cell lysate is kept in -20°C or -80°C. You should process the samples all the way to isopropanol precipitation/ethanol washes (even though you may leave the isopropanol-filled tubes at -80°C for weeks/months).



Assessment of RNA Quality and Integrity Is Essential for Downstream Analyses (esp. RT-PCR and qPCR). (A) A total RNA sample was degraded for varying times at -20C and the resulting samples were analyzed on the Agilent 2100 Bioanalyzer System using the Eukaryote Total RNA Nano assay. A shift towards shorter fragment sizes can be observed with progressing degradation. (B) Sample electropherograms used to train the RNA Integrity Number (RIN) software. Samples range from intact (RIN 10), to degraded (RIN 2). (C) Representative “good” and “bad” RNA samples.